



HA Cartridges FAQ

📍 Address: No. 98, Technology Sixth Road, High-tech Zone, Zhuhai City, 519085, Guangdong, China

☎ Tel: +86 (756) 3689708

✉ E-mail: customerservice@jafron.com

🌐 Website: www.jafron.com



Basic FAQ Related to Jafron Hemoperfusion

1. How do the structure of Jafron HA resins look like?

- The adsorbent materials of HA cartridges are neutro-macroporous resin made of double cross-linked styrene-divinylbenzene copolymers.
- The adsorption principles of the resins include the molecular sieve, van der waals, and the hydrophobic characteristics of the resins.
- Jafron's porosity control technology adjusts the resin pore size distribution and adsorption range.
- Collodion coats the resins to mitigate the possible adverse effect of direct contact with blood.

2. What are the indications for HA cartridges?

According to clinical practices and studies, disposable hemoperfusion cartridge is indicated to remove the following substances^[1-11]:

- Inflammatory mediators and cytokines such as IL-1, IL-6, IL-8, IL-10, TNF- α .
- Overdosed drugs and poisons such as organophosphorus, paraquat, carbamazepine.
- Accumulated β 2-MG, PTH, leptin and protein-bound toxins in end stage renal disease hemodialysis-related complication.
- Excessive triglyceride and cholesterol in hyperlipidaemia severe acute pancreatitis.
- Other substances: bilirubin, myoglobin.

3. Could you illustrate the safety and biocompatibility of Jafron cartridges?

The biocompatibility of the cartridges has been evaluated both in vitro and in vivo. In vitro, CT imaging showed an excellent distribution of the flow inside the HA cartridges without channeling phenomena^[10]. The cytotoxicity test demonstrated that HA cartridges carried out an optimal level of biocompatibility, and their use was not associated with adverse reactions^[11]. Clinically, the publications did not show significant adverse reactions related to Jafron hemoperfusion^[1,2,8,12,13].

4. What is the platelet loss rate during hemoperfusion with Jafron?

According to the previous studies, Jafron hemoperfusion did not show a severe reduction in platelet levels^[14-15]; the loss of the platelets is low, usually is less than 10%, and usually recovered within 24 hours.

5. Do hemoglobin and albumin adsorbed by the cartridges?

The previous studies did not show significant decrease of hemoglobin and albumin after the hemoperfusion treatment^[16,17].

- [1] Huang, Zhao, et al. "Effect on extrapulmonary sepsis - induced acute lung injury by hemoperfusion with neutral microporous resin column." *Therapeutic Apheresis and Dialysis* 17.4 (2013): 454-461.
- [2] Sun, Shiren, et al. "High-volume hemofiltration plus hemoperfusion for hyperlipidemic severe acute pancreatitis: a controlled pilot study." *Annals of Saudi medicine* 35.5 (2015): 352-358.
- [3] Xu, Xuefeng, et al. "Effect of HA330 resin-directed hemoabsorption on a porcine acute respiratory distress syndrome model." *Annals of intensive care* 7.1 (2017): 1-17.
- [4] Li, Li, et al. "Hemoperfusion plus continuous veno-venous hemofiltration in the treatment of patients with multiple organ failure after wasp stings." *The International journal of artificial organs* 43.3 (2020): 143-149.
- [5] Li, An, et al. "Early stage blood purification for paraquat poisoning: a multicenter retrospective study." *Blood purification* 42.2 (2016): 93-99.
- [6] Jiang, Shu-zhi, et al. "Clinical efficacy of intravenous infusion of atropine with micropump in combination with hemoperfusion on organophosphorus poisoning." *Saudi journal of biological sciences* 26.8 (2019): 2018-2021.
- [7] Yang, Xiangming, et al. "Early hemoperfusion for emergency treatment of carbamazepine poisoning." *The American journal of emergency medicine* 36.6 (2018): 926-930.
- [8] Huu, Dung Nguyen, et al. "A Combination of Hemodialysis with Hemoperfusion Helped to Reduce the Cardiovascular-Related Mortality Rate after a 3-Year Follow-Up: A Pilot Study in Vietnam." *Blood purification* 50.1 (2021): 65-72.
- [9] Chen, Shun-Jie, et al. "Combination of maintenance hemodialysis with hemoperfusion: a safe and effective model of artificial kidney." *The International journal of artificial organs* 34.4 (2011): 339-347.
- [10] Lorenzin, A., et al., Fluid Dynamics Analysis by CT Imaging Technique of New Sorbent Cartridges for Extracorporeal Therapies. 2019. 48(1): p.18-24.
- [11] Montin, D.P., et al., Biocompatibility and cytotoxic evaluation of new sorbent cartridges for blood hemoperfusion. 2018. 46(3): p. 187-195.
- [12] Qing-hua, W., and W. Feng. "Effect of Coupled Plasma Filtration Adsorption on Inflammatory Mediators and Liver Function of Patients with Severe Acute Pancreatitis 2019; 2: 111." *Research Article*.
- [13] Kidney." *The International journal of artificial organs* 34.4 (2011): 339-347.
- [14] Yuan, Hai, et al. "Efficacy of two combinations of blood purification techniques for the treatment of multiple organ failure induced by wasp stings." *Blood purification* 42.1 (2016): 49-55.
- [15] Kaçar CK, Uzundere O, Kandemir D, Yekta A. Efficacy of HA330 Hemoperfusion Adsorbent in Patients Followed in the Intensive Care Unit for Septic Shock and Acute Kidney Injury and Treated with Continuous Venovenous Hemodiafiltration as Renal Replacement Therapy. *Blood Purif.* 2020;49(4):448-456. doi: 10.1159/000505565. Epub 2020 Jan 28. PMID: 31991412.
- [16] Yuan, Hai, et al. "Efficacy of two combinations of blood purification techniques for the treatment of multiple organ failure induced by wasp stings." *Blood Purification* 42.1 (2016): 49-55.
- [17] Chen, Shun-Jie, et al. "Combination of maintenance hemodialysis with hemoperfusion: a safe and effective model of artificial kidney." *The International journal of artificial organs* 34.4 (2011): 339-347.

Contraindications, warnings and precautions of the products, please refer to Instruction For Use.

HA130 Related FAQ

1. What are the main clinical applications of HA130 hemoperfusion?

MHD patients have the following clinical manifestations, they recommended to the HP treatment^[1]:

- Severe uremic pruritus
- Severe uremic-related sleep disorders
- Protein energy consumption
- Microinflammatory status
- Severe secondary hyperparathyroidism
- Severe hyper β 2-MG
- Refractory hypertension
- Restless Leg Syndrome (RLS)
- Uremic peripheral neuropathy

2. What is the treatment scheme for Ha130?

The treatment scheme of HA130 could be as follow^[2,3]:

- Intensive program: 4 times/month
- Maintenance program: 1-2 times/month
- Treatment duration: 4 hours combined with HD or HDF
- The treatment scheme could be adjusted according to the patient's condition.

3. In some cases, why is the PTH still higher or did not reduce after the HA130 treatment?

HP treatment is expected to improve the symptoms and reduce the middle uremic toxins and protein-bound uremic toxins; try to analyze the factors which may affect the results:

- Patient's diet: A high phosphate diet and phosphate levels may increase PTH secretion^[4].
- Serum low vitamin D and calcium may increase the phosphate levels and increase PTH secretion^[4].

- Secondary hyperparathyroidism: Associated with high levels of PTH^[5].
- Note: Compare the HA130 group to the non-HA130 group to better outline the effect of HA130.

4. Is there any impact of HA130 on quality of life?

Combining HA130 and HD could improve body pain, general health, vitality, emotional role, mental health, and total QoL score^[2].

5. Is there an impact of HA130 on survival?

The Systematic Review and Meta-Analysis on the survival outcomes for the patients with End-Stage Renal Disease showed that the 1-year overall survival (OS) rate, the 2-year OS rate, and the 5-year OS rate of patients with ESRD treated with HP + HD were better than those treated with HD^[6].

6. Is there a cost-effectiveness study on hemoperfusion treatment?

The cost-effectiveness study with 1,407 patients recruited to the HD/HP trial from 30 clinical centers with two years' follow-up period demonstrated that^[7]:

- In the base case analysis, compared with HD alone, HD + HP results in 2.87 LYs saved, 1.32 QALY gains and an additional cost of RMB 33,340 per patient. The probabilistic ICER of HD + HP is USD 25,251 per QALY, which is lower than the USD 30,778 willingness-to-pay threshold. Therefore, HD + HP is considered to be cost-effective.
- HP+HD reduces the incidence of severe CVD events and subsequent CVD deaths.

[1] Gengru Jiang, "Shanghai Expert Consensus of Hemoperfusion Therapy Application in Maintenance Hemodialysis Patients", 39th Vicenza Course on AKI & CRRT, 26-29 October 2021.

[2] Chen, Shun-Jie, et al. "Combination of maintenance hemodialysis with hemoperfusion: a safe and effective model of artificial kidney." *The International journal of artificial organs* 34.4 (2011): 339-347.

[3] Huu, Dung Nguyen, et al. "A Combination of Hemodialysis with Hemoperfusion Helped to Reduce the Cardiovascular-Related Mortality Rate after a 3-Year Follow-Up: A Pilot Study in Vietnam." *Blood purification* 50.1 (2021): 65-72.

[4] Takeda, Eiji, et al. "Increasing dietary phosphorus intake from food additives: potential for negative impact on bone health." *Advances in nutrition* 5.1 (2014): 92-97.

[5] Yuen, Noah K., Shubha Ananthakrishnan, and Michael J. Campbell. "Hyperparathyroidism of renal disease." *The Permanente Journal* 20.3 (2016).

[6] Cheng, Wendi, et al. "Survival Outcomes of Hemoperfusion and Hemodialysis versus Hemodialysis in Patients with End-Stage Renal Disease: A Systematic Review and Meta-Analysis." *Blood Purification* (2021): 1-13.

[7] Wang, Haiyin, et al. "Cost-effectiveness analysis of hemodialysis plus hemoperfusion versus hemodialysis alone in adult patients with end-stage renal disease in China." *Annals of Translational Medicine* 9.14 (2021).

Contraindications, warnings and precautions of the products, please refer to Instruction For Use.

HA230 Related FAQ

1. What are the adsorption principles and the clinical applications of the HA230 hemoperfusion cartridge?

It mainly removes poisons and overdose drugs with a molecular weight of 500Da~10KDa or hydrophobic or high protein-binding properties^[1].

The clinical applications could be divided into poisons and drug overdose as listed below^[2-12] :

- Parquat poisoning, Acute severe organophosphate poisoning (ASOP), 2,4-Dinitrophenol, Snake venom, Mushroom poisoning, Abrin, Paroxetine, Thallium, Aconitine, etc.
- Methotrexate (MTX), Ticagrelor, Carbamazepine, Paroxetine, Amitriptyline, Digoxin, Acetaminophen, Imipramine, Amlodipine (CCBs), etc.

2. What is the recommended treatment scheme for HA230?

- The treatment should be initiated as early as possible (4-6 h from exposure to poison or drug, the therapy after 12 hours of the exposure is less effective).
- The recommended scheme is 3-5 cartridges for each patient, with a suggested duration of 3-4 hours for each cartridge and no more than 12 hours.
- Hemoperfusion should be carried out more frequently at the beginning (1-3 cartridges/day).
- The treatment frequency and duration could be adjusted according to the type of poison/drug and the patient's condition.

3. What are the expected benefits of using HA230 hemoperfusion?

According to the previous publications, hemoperfusion HA230 therapy could benefit on^[2,11,12]:

- Accelerate the elimination of poison and overdose drugs.
- Stabilize the hemodynamic status.

- Alleviate complications such as consciousness disorder, respiratory depression, coma, etc.
- Shorten hospitalization.
- Improve the survival.

4. What are the treatment modes for HA230 hemoperfusion cartridge?

HA230 could stand alone as single hemoperfusion or combined with HD or CRRT.

- [1] Ankawi, Ghada, et al. "A new series of sorbent devices for multiple clinical purposes: current evidence and future directions." *Blood purification* 47.1-3 (2019): 94-100.
- [2] L.BO, et al. "Therapeutic efficacies of different hemoperfusion frequencies in patients with organophosphate poisoning". *European Review for Medical and Pharmacological Sciences* 2014; 18:3521-3523.
- [3] Hui Dong, Bachelor, Yi-Bing Weng, MD, Gen-Shen Zhen, Bachelor, Feng-Jie Li, MM, Ai-Chun Jin, Bachelor, Jie Liu, MM, et al. "Clinical emergency treatment of 68 critical patients with severe organophosphorus poisoning and prognosis analysis after rescue". *Medicine* (2017) 96:25(e7237)
- [4] Deven Juneja, Omender Singh, Alka Bhasin, Manish Gupta, Sanjay Saxena, Archana Chaturvedi, et al. "Severe Suicidal Digoxin Toxicity Managed with Resin Hemoperfusion: A Case Report". *Indian Journal of Critical Care Medicine*. 2011
- [5] Xue-hong ZHAO, Jiu-kun JIANG, Yuan-qiang LU, et al. "Evaluation of efficacy of resin hemoperfusion in patients with acute 2,4-dinitrophenol poisoning by dynamic monitoring of plasma toxin concentration". *Journal of Zhejiang University-SCIENCE B (Biomedicine & Biotechnology)* ISSN 1673-1581 (Print); ISSN 1862-1783 (Online)
- [6] Jiliang Huang, MD, Wenbin Zhang, MD, Xin Li, MD, Shufen Feng, MD, et al. "Acute abrin poisoning treated with continuous renal replacement therapy and hemoperfusion successfully". *Medicine* (2017) 96:27(e7423).
- [7] Junxiu Zhao, Xiaobo Peng, Chunyan Bai, Lili Bai, et al. "Efficacy analysis of prussian blue or its combination with hemoperfusion in the treatment of acute thallium". *Chin Crit Care Med*, July 2018, Vol.30, No.7.
- [8] Yanling Chen, Liang Li, Feng Sheng, Wei Xiao, et al. "Clinical experience of using continuous renal replacement therapy combined with hemoperfusion successfully". *Chin J TCM WM Crit Care*, November 2018, Vol.25, No.6.
- [9] Siyuan Jiang, Wei Ke, Fu Zhang, et al. "Clinical study on blood perfusion combined with naloxone for rescuing patients with hypnotic poisoning". *China Pharmaceuticals*, June, 5, 2019, Vol.28, No.11.
- [10] Meihui Jiang, Weimin Zheng, Rui Fu, Yuhui, Huang, et al. "Clinical analysis on hemoperfusion combined with hemodialysis in the treatment of children with severe mushroom poisoning". *Jiangxi Medical Journal*, December 2017, Vol 52, No12.
- [11] Xiangming Yang, Shiyu Xin, Yajie Zhang, Tiegang Li, et al. "Early hemoperfusion for emergency treatment of carbamazepine poisoning". *ELSEVIER YAJEM-57051*: No of Pages 5.
- [12] An Li, Wenxiong Li, Fengtong Hao, Haishi Wang, et al. " Early Stage Blood Purification for Paraquat Poisoning: A Multicenter Retrospective Study." *Blood Purif* 2016;42:93–99; DOI: 10.1159/000445991.

Contraindications, warnings and precautions of the products, please refer to Instruction For Use.

HA380 Related FAQ

1. What are the main clinical applications of Ha380?

- Sepsis/Septic shock^[1,2]
- COVID-19^[3,4]
- Severe acute pancreatitis/ hyperlipidemia induced pancreatitis^[5]
- Trauma/Burn^[6]
- Cardiopulmonary bypass surgery^[7]

2. When to initiate the HA380 hemoadsorption treatment?

The timing of the application of hemoperfusion for critical illness is still under investigation. The inclusion criteria from previous publications and expert experiences may differ in different countries; however, inflammatory mediators of critically ill patients may release in the early stage of the disease. The early use could have more benefits and control the disease's development.

At the 2021 International Critical Illness Conference (ISICEM), Professor Claudio Ronco and others stated that:

- Start early to remove the pro-inflammatory cytokines in the first 24 hours of admission to the ICU and before the organ dysfunction.
- Start immediately for severe burn patients and long period cardiac bypass CPB patients.
- For other types of patients, especially during the coronavirus, do not wait for the indication of CRRT to perform hemoperfusion.
- A recent publication demonstrated that the early start of hemoperfusion could be more effective and significantly reduce the mortality rate among COVID-19 patients with critical diseases^[3].

3. What is the recommended treatment scheme for HA380?

1~2 cartridges/day for three consecutive days, the suggested duration of each cartridge is 6 hours, no more than 12 hours. The physician could adjust the treatment scheme according to the patient's condition.

4. What are the expected benefits of using the HA380 hemoperfusion cartridge?

HA380 hemoperfusion may benefit the patients by^[2-8]:

- Reduce inflammatory cytokines and other target substances such as IL-6, IL-8, IL-1 and TNF- α .
- Reduce the triglyceride and cholesterol levels of severe acute pancreatitis.
- Improve hemodynamics and reduce the vasopressor requirement.
- Improve organ function, oxygenation index, SOFA score, APACHE II score.
- Reduce mechanical ventilation, shorten ICU stay, and improve the survival.

5. Could HA380 be used for traumatic patients?

Hemoperfusion could reduce myoglobin, creatinine kinase and stabilize the hemodynamics in indicated patients^[9].

6. Could HA380 be used for severe acute pancreatitis (SAP)?

Hemoperfusion could reduce the cytokines and lipids (Triglyceride & cholesterol). It could reduce the severity of the disease and stabilize the hemodynamics^[10].

7. Could HA380 adsorb the Ticagrelor?

The lab results showed that hemoperfusion could remove the Ticagrelor, and the removal rate could reach 80% in 6h^[11].

[1] Huang, Zhao, et al. "Effect on Extrapulmonary Sepsis - Induced Acute Lung Injury by Hemoperfusion with Neutral Microporous Resin Column." Therapeutic Apheresis and Dialysis 17.4 (2013): 454-461.

[2] Arslan, Baris, et al. "A single-center experience with resin adsorption hemoperfusion combined with continuous veno-venous hemofiltration for septic shock patients." Med Science 8 (2019): 390-4.

[3] Mikaeili, Haleh, et al. "The early start of hemoperfusion decreases the mortality rate among severe COVID - 19 patients: A preliminary study." Hemodialysis International.

[4] Dastan, Farzaneh, et al. "Continues renal replacement therapy (CRRT) with disposable hemoperfusion cartridge: a promising option for severe COVID-19." Journal of global antimicrobial resistance 21 (2020): 340-341.

[5] Sun, Shiren, et al. "High-volume hemofiltration plus hemoperfusion for hyperlipidemic severe acute pancreatitis: a controlled pilot study." Annals of Saudi medicine 35.5 (2015): 352-358.

[6] Li, Li, et al. "Hemoperfusion plus continuous veno-venous hemofiltration in the treatment of patients with multiple organ failure after wasp stings." The International Journal of Artificial Organs 43.3 (2020): 143-149.

[7] Zijian He, et al. Efficacy of Resin Hemoperfusion Cartridge on Inflammatory Responses during Adult Cardiopulmonary Bypass. Blood Purif. 2021 Jun 9;1-7.

[8] Chavez, Joselito R., et al. "A case of leptospirosis with acute respiratory failure and acute kidney injury treated with simultaneous extracorporeal membrane oxygenation and haemoperfusion." BMJ Case Reports CP 12.5 (2019): e229582.

[9] Yuan, Hai, et al. "Efficacy of two combinations of blood purification techniques for the treatment of multiple organ failure induced by wasp stings." Blood purification 42.1 (2016): 49-55.

[10] Sun, Shiren, et al. "High-volume hemofiltration plus hemoperfusion for hyperlipidemic severe acute pancreatitis: a controlled pilot study." Annals of Saudi medicine 35.5 (2015): 352-358. plus hemoperfusion versus hemodialysis alone in adult patients with end-stage renal disease in China." Annals of Translational Medicine 9.14 (2021).

[11] Du Hongyan et al., "Ticagrelor Removal by Hemoadsorption From Human Blood", 39th Vicenza Course on AKI & CRRT, 26-29 October 2021.

Contraindications, warnings and precautions of the products, please refer to Instruction For Use.

DPMAS Related FAQ

1. Why do we need bilirubin adsorption?

High bilirubin levels could lead to mitochondrial dysfunction and capillary leak syndrome; cholemic nephrosis is also associated with elevated bilirubin in ACLF. Bile acid could lead to hepatotoxicity and nephrotoxicity. Systemic inflammation and bile acids are critical drivers for liver failure in ACLF. High ammonia is supposed to be associated with brain swelling. DPMAS is not only to adsorb the bilirubin but also to adsorb bile acid, ammonia, and cytokines; this may benefit on:

- Decrease bilirubin, bile acid, cytokines, and ammonia^[1-5].
- Alleviate liver injury^[1-6].
- Reduce hepatic encephalopathy grading^[6].
- Increase liver transplant waiting time^[2,3].

2. What are the treatment modes for liver cartridges?

Jafron hemoadsorption therapy for the liver disease could be performed in different modes:

- HA330-II: Direct blood adsorption as HP, HP+HD or HD+CRRT mode.
- BS3330: Plasma adsorption (PA) mode (Plasma separator + BS330 cartridge).
- DPMAS: Plasma adsorption (PA) mode (Plasma separator+BS330 cartridge+HA330-II cartridge).

3. When do we use HA330-II, BS330 and DPMAS?

HA330-II	BS330	DPMAS
Hyper-inflammation + Hyperbilirubinemia	Hyperbilirubinemia + Hyper-bileacidemia	Hyperbilirubinemia + Hyper-inflammation
Total Bilirubin < 13mg/dl; High PCT, CRP, or IL-6	Total Bilirubin > 222.4 umol/L (13mg/dl).	Total Bilirubin ≥ 13mg/dl; High PCT,CRP, or IL-6

4. Why do we need a plasma separator when using BS330?

The direct contact of blood to BS330 resin may cause severe adverse reactions; the reasons are:

- The BS330 resins are anion exchange resins; direct exposure to the blood may cause coagulation.
- BS330 cartridge has finer resins than HA cartridges; whole blood adsorption may lead to cartridge blockage.
- Filtration fabric at the outlet end of the BS330 column may result in clotting if direct blood adsorption.

5. Why sometimes there are white or translucent jelly-like substances appear on BS330 during the process of plasma adsorption?

- It is a sign of coagulation problem resulting from inadequate anticoagulation. There is no blood cell, so it will not form a clot. The visible white or translucent jelly-like material is activated fibrinogen and precipitation phenomena.
- Mostly, it is caused by insufficient anticoagulation, or the patient is in a hypercoagulable state. Arterial pressure, venous pressure, and TMP should be monitored carefully during the treatment. The adjustment of anticoagulant should be considered.

6. When combining PE+DPMAS, which treatment should start first?

- PTA<30%, recommend conducting PE then DPMAS.
- PTA>30%, could start from either PE or DPMAS.
- PE is recommended to be admitted first for patients with PTA<30%, because it could supply the patients with the coagulation factors and improve the coagulation function.

[1] Molecular absorption system in patient with decompensated liver failure." *nephrology*. Vol. 25. 111 river st, hoboken 07030-5774, nj usa:1] Yao, Jia, et al. "Therapeutic effect of double plasma molecular adsorption system and sequential half - dose plasma exchange in patients with HBV - related acute - on - chronic liver failure." *Journal of clinical apheresis* 34.4 (2019): 392-398.

[2] Li, Peng, et al. "A non-bioartificial liver support system combined with transplantation in HBV-related acute-on-chronic liver failure." *Scientific reports* 11.1 (2021): 1-9.

[3] Wan, Yue - Meng, et al. "Therapeutic plasma exchange versus double plasma molecular absorption system in hepatitis B virus - infected acute - on - chronic liver failure treated by entecavir: A prospective study." *Journal of clinical apheresis* 32.6 (2017): 453-461.

[4] Guo, Xiju, et al. "Comparison of plasma exchange, double plasma molecular adsorption system, and their combination in treating acute-on-chronic liver failure." *Journal of International Medical Research* 48.6 (2020): 0300060520932053.

[5] Sharma, Dhruva, et al. "Hepatitis A virus-induced severe hemolysis complicated by severe glucose-6-phosphate dehydrogenase deficiency." *Indian journal of critical care medicine: peer-reviewed, official publication of Indian Society of Critical Care Medicine* 22.9 (2018): 670.

[6] Sirivongrangsorn, Phatadon, et al. "immunomodulatory effect of double plasma wiley, 2020.

Contraindications, warnings and precautions of the products, please refer to Instruction For Use.

HA Cartridges Operation Related FAQ

1. Why should we heparinize the HA cartridges?

- Heparin molecular weight: 15kDa.
- Heparin is expected to be adsorbed if used for systemic anticoagulation.
- Jafron recommends heparinization by injecting heparin into the cartridge 30 minutes before priming to prevent the possible adsorption of heparin and prevent coagulation problem, ensuring more safety procedures.
- Citrate anticoagulation is recommended for patients who do not tolerate heparinization.

2. Should we heparinize the cartridge when using citrate anticoagulation?

If the patient will go through citrate anticoagulation, you could prime the HA cartridge with 3000mL normal saline without heparin.

3. What is the blood flow rate and the duration of HP treatment?

The recommended blood flow rate and treatment duration are as follow:

HP	Blood flow rate: increase gradually from 100~150mL/min to 200~250mL/min. Duration: 2~2.5 hours.
HP+HD/HDF	Blood flow rate: usually less than 320mL/min. Duration: 4 hours.
HP+CRRT	Blood flow rate: usually 150~250mL/min. Duration: can be up to 12 hours.
HP+CPB	Blood flow rate: can be up to 700mL/min. Duration: usually less than 2.5 hours.
Blood flow rate: 100 ~700mL/min. Treatment duration: 2~12 hours.	

4. What kind of machines do we need to perform the hemoperfusion treatment with Jafron cartridges?

Hemoperfusion using HA cartridges is compatible with all common hemodialysis and CRRT machines.

5. How do you anticoagulate the patients during the hemoperfusion treatment?

The suggested anticoagulant type and dosage are as follow:

Heparin

Initial recommended dosage: 62.5~125U/kg.

Additional dosage:1250~2500U/hour.

Stop adding heparin 30 minutes before the end of treatment.

Low molecular weight heparin

It might be used for anticoagulation according to treatment requirements.

Citrate

When using regional anticoagulation, the concentration of citrate and calcium ions should be closely monitored during the entire treatment.

***Note:**

Due to the individual differences, above anticoagulant dosage should depend on patient's condition.

Other anticoagulation agents may be used according to the patient's condition.

Contraindications, warnings and precautions of the products, please refer to Instruction For Use.



- Established in 1989, IPO on ChiNext in 2016
- Widely used in more than 80 countries
- Applied in over 8000 influential hospitals
- Annual application over 5 million pcs



*Contraindications, Warnings and Precautions refer to Instructions For Use

JAFRON BIOMEDICAL CO., LTD.

Address: No. 98, Technology Sixth Road, High-tech Zone, Zhuhai City, 519085, Guangdong, China.

Tel: +86 (756) 3689708

E-mail: customerservice@jafron.com

Website: www.jafron.com

(For Internal Use)